



DESIGN / ENGINEER / MANUFACTURE

## Barrier-First - Product Specification Community Face Covering - Made in Scotland

### Innovative Safety

- Highly effective melt blown polypropylene filtration layer with spun bond polypropylene outer layers
- UV sterilized
- Automatically hygienically wrapped
- Universal Fit
- Single Use

| Finished Covering  |   |
|--|---|
| Unit Weight  | Target: 3 g/m <sup>2</sup>  |
| Colour   | White   |
| Dimensions<br>(Standard Type 1/II/IIR medical mask Dimensions) | (L) 175mm +/- 5mm<br>(W): 95mm +/- 5mm<br>Ear Loop 180mm +/- 10mm |

|  |
|--|
| <b>UK Sourced Materials - 3 Layer Non-Woven Construction</b> |
| Spun bond polypropylene inner and outer layer                |
| Melt Blown Polypropylene Filter Layer                        |
| Latex Free nylon Spandex Ear Loop                            |
| Integrated PE/Zinc Galvanised Nose Wire                      |

|  |   |
|--|---|
| Conforms to AFNOR S76 001  | Inhalation @ 30l/min<br>Inhalation @ 95 l/min<br>Exhale @ 160 l/min<br>Filtration (NaCl Penetration) >70% |
| Manufactured under stringent quality controls and in compliance with ISO 9001 regulations guaranteeing quality, performance and ensuring full material and product traceability. |   |

Transcal's Community Face Coverings are to be used in conjunction with the latest government Covid-19 guidelines and in line with social distancing regulations, more information can be found on the government website ([www.gov.uk](http://www.gov.uk)).



T: +44 (0)1506 440 111  
F: +44 (0)1506 442 333  
E: [info@transcal.co.uk](mailto:info@transcal.co.uk)  
W: [www.transcal.co.uk](http://www.transcal.co.uk)

**TRANSCAL**

TRANSCAL LIMITED  
VAT NO. GB 446 4182 46  
REG IN SCOTLAND NO. 96703



DESIGN / ENGINEER / MANUFACTURE



BARRIER FIRST

### Guidance for Use

Face coverings should be fitted correctly ensuring your nose and chin are covered allowing you to breathe comfortably.

- Wash your hands thoroughly with soap and water for 20 seconds or use hand sanitizer before putting on a face covering.
- The face covering should fit comfortably but securely against the side of the face.
- Avoid touching the part of the face covering in contact with your mouth and nose.
- Secure behind the ears with the elastic ear loops.
- Before removing wash your hands with soap and water for 20 seconds or use hand sanitizer.
- Only handle the ear loops.
- Dispose of the face covering in a safe and responsible manner, not recyclable.
- Once disposed of wash your hands with soap and water for 20 seconds or use hand sanitizer.

### Further Information

Change covering if it becomes damp or if you have touched it.

Avoid taking it off and putting it back on a lot in quick succession (for example, when leaving and entering shops on a high street).

Do not give your face covering to anyone else to use.

More guidance can be obtained from [www.gov.uk](http://www.gov.uk)

Date: 27<sup>th</sup> August 2020



Certificate No. 522  
ISO 9001:2008



UK.21G.2653

T: +44 (0)1506 440 111  
F: +44 (0)1506 442 333  
E: [info@transcal.co.uk](mailto:info@transcal.co.uk)  
W: [www.transcal.co.uk](http://www.transcal.co.uk)

**TRANSCAL**

TRANSCAL LIMITED  
VAT NO. GB 446 4182 46  
REG IN SCOTLAND NO. 96703